

Glyt2 Rabbit mAb

ABZ1736


Product Overview

Name	Glyt2 Rabbit mAb
Accession(Primary)	Q9Y345
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9Y345)

Gene Name SLC6A5

About Gene Symbol



Western blot analysis of Glyt2 in Jurkat, C6 lysates using GlyT2 antibody.

