

AP3M1 Rabbit mAb


ABZ1738

Product Overview

Name	AP3M1 Rabbit mAb
Accession(Primary)	Q9Y2T2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q9Y2T2)

Gene Name	AP3M1
About Gene Symbol	AP3M1 (Adaptor Related Protein Complex 3 Subunit Mu 1) is a signaling adaptor protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to eluc



Western blot analysis of AP3M1 in K562, C6, 3T3, HeLa lysates using AP3M1 antibody.