

SSX2IP Rabbit mAb


ABZ1740

Product Overview



Name	SSX2IP Rabbit mAb
Accession(Primary)	Q9Y2D8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; IP-1:100-1:200; ICC/IF-1:50-1:100

Target information (Q9Y2D8)


Gene Name	SSX2IP
About Gene Symbol	




Immunocytochemistry analysis of SSX2IP (green) in K562 using SSX2IP antibody, and DAPI (blue).



Western blot analysis of SSX2IP in rat kidney lysates using SSX2IP antibody.



Western blot analysis of SSX2IP in Jurkat, HeLa lysates using SSX2IP antibody



Immunohistochemistry analysis of paraffin-embedded Human Brain using SSX2IP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.