

Phospholipase A2 Activator Protein Rabbit mAb


ABZ1742

Product Overview


Name	Phospholipase A2 Activator Protein Rabbit mAb
Accession(Primary)	Q9Y263
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ICC/IF-1:50-1:100

Target information (Q9Y263)

Gene Name	PLAA
About Gene Symbol	PLAA (Phospholipase A2 Activating Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Phospho-lipase A2 Activator Protein (green) in Hela using Phospho-lipase A2 Activator Protein antibody, and DAPI (blue).



Western blot analysis of Phospho-lipase A2 activator protein (PLAP) in 3T3, Hela lysates using Phospho-lipase A2 activator protein (PLAP) antibody.