

AKT Rabbit mAb


ABZ1743

Product Overview


Name	AKT Rabbit mAb
Accession(Primary)	Q9Y243
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:500-1:1000

Target information (Q9Y243)


Gene Name	AKT3
About Gene Symbol	AKT3 (AKT Serine/Threonine Kinase 3) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of AKT in 3T3, rat heart lysates using AKT antibody.



Western blot analysis of AKT1/2/3 in LnCap lysates using AKT antibody.



Immunocytochemistry analysis of
AKT1/2/3(green) in Hela using AKT antibody, and
DAPI(blue)