


GNE Rabbit mAb

ABZ1744

Product Overview

Name	GNE Rabbit mAb
Accession(Primary)	Q9Y223
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20; IF-1:50-1:200

Target information (Q9Y223)

Gene Name GNE**About Gene Symbol**

Western blot analysis of GNE in K562, HeLa lysates using GNE antibody.

