

STAT2 Rabbit mAb

ABZ1745


Product Overview

Name	STAT2 Rabbit mAb
Accession(Primary)	Q9WVL2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 105 kDa; Observed MW: 113 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (Q9WVL2)

Gene Name Stat2

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded mouse lung using STAT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of STAT2 in HL-60, U251, U87-MG lysates using STAT2 antibody.