

CaMKII alpha Rabbit mAb


ABZ1746

Product Overview


Name	CaMKII alpha Rabbit mAb
Accession(Primary)	Q9UQM7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC:1:200-1:500; ICC/IF-1:20-1:100; FC-1:20-1:100

Target information (Q9UQM7)

Gene Name	CAMK2A
About Gene Symbol	



Western blot analysis of CaMKII alpha in 3T3 lysates using CaMKII alpha antibody.



Western blot analysis of CaMKII alpha in rat Brain, HeLa lysates using CaMKII alpha antibody.