


xCT Rabbit mAb

ABZ1747

Product Overview

Name	xCT Rabbit mAb
Accession(Primary)	Q9UPY5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:4000; IF-1:50-1:200; IP-1:50

Target information (Q9UPY5)

Gene Name SLC7A11**About Gene Symbol**

Western blot analysis of xCT in K562 lysates using xCT antibody.

