

SATB2 Rabbit mAb


ABZ1748

Product Overview

Name	SATB2 Rabbit mAb
Accession(Primary)	Q9UPW6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:200-1:1000


Target information (Q9UPW6)

Gene Name	SATB2
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.


Western blot analysis of SATB2 in HT-1080, Rat brain, Mouse brain lysates using SATB2 antibody.




Western blot analysis of SATB2 in HT1080 Rat brain lysates using SATB2 antibody.



Western blot analysis of SATB2 in Saos2 lysates using SATB2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human stomach cancer using SATB2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of SATB2 in SH-SY5Y using SATB2 antibody.