


USP22 Rabbit mAb

ABZ1749


Product Overview

Name	USP22 Rabbit mAb
Accession(Primary)	Q9UPT9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20; IHC-1:50-1:100

Target information (Q9UPT9)

Gene Name USP22**About Gene Symbol**

Western blot analysis of USP22 in F9 lysates using USP22 antibody.



Western blot analysis of USP22 in A549, HL-60, U2OS, C6 lysates using USP22 antibody.