

## STUB1 Rabbit mAb


ABZ1752

### Product Overview


|                             |                                                      |
|-----------------------------|------------------------------------------------------|
| <b>Name</b>                 | STUB1 Rabbit mAb                                     |
| <b>Accession(Primary)</b>   | Q9UNE7                                               |
| <b>Host Species</b>         | Rabbit                                               |
| <b>Clonality</b>            | Monoclonal                                           |
| <b>Isotype</b>              | IgG                                                  |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                    |
| <b>Calculated MW</b>        | Calculated MW: 35 kDa; Observed MW: 35 kDa           |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; ICC/IF-1:50-1:200;<br>IP-1:20-1:50 |

### Target information (Q9UNE7)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | STUB1                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | STUB1 (STIP1 Homology And U-Box Containing Protein 1) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and |



Western blot analysis of STUB1 in K562, rat Brain, C6, 3T3, HeLa lysates using STUB1 antibody.



Immunocytochemistry analysis of STUB1 (green) in HeLa using STUB1 antibody, and DAPI (blue)