

VPS4A Rabbit mAb


ABZ1753

Product Overview


Name	VPS4A Rabbit mAb
Accession(Primary)	Q9UN37
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100

Target information (Q9UN37)


Gene Name	VPS4A
About Gene Symbol	



Western blot analysis of VPS4a in K562, rat Brain, C6, 3T3 lysates using VPS4A antibody.



Immunocytochemistry analysis of VPS4a(green) in HeLa using VPS4a antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using VPS4a antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.