

Bok Rabbit mAb

ABZ1755


Product Overview

Name	Bok Rabbit mAb
Accession(Primary)	Q9UMX3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:20-1:50

Target information (Q9UMX3)

Gene Name BOK

About Gene Symbol



Western blot analysis of Bok in rat Brain, C6, HeLa lysates using Bok antibody.

