

TPX2 Rabbit mAb


ABZ1756

Product Overview


Name	TPX2 Rabbit mAb
Accession(Primary)	Q9ULW0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 86 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q9ULW0)


Gene Name	TPX2
About Gene Symbol	TPX2 (TPX2 Microtubule Nucleation Factor) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targ



Western blot analysis of TPX2 in HeLa, A549, HL-60 lysates using TPX2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TPX2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TPX2 (green) in Hela using TPX2 antibody, and DAPI(blue)