

RAB22A Rabbit mAb


ABZ1757

Product Overview

Name	RAB22A Rabbit mAb
Accession(Primary)	Q9UL26
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; IP-1:10-1:100

Target information (Q9UL26)

Gene Name	RAB22A
About Gene Symbol	



Western blot analysis of RAB22A in rat Brain, C6, 3T3, HeLa lysates using Rab22A antibody.