

NUP50 Rabbit mAb


ABZ1759

Product Overview


Name	NUP50 Rabbit mAb
Accession(Primary)	Q9UKX7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (Q9UKX7)

Gene Name	NUP50
About Gene Symbol	



Immunocytochemistry analysis of NUP50 (green) in Jurkat using NUP50 antibody, and DAPI (blue).



Western blot analysis of NUP50 in HeLa, A549, HL-60, U251 lysates using NUP50 antibody.