

PACSIN3 Rabbit mAb


ABZ1762

Product Overview


Name	PACSIN3 Rabbit mAb
Accession(Primary)	Q9UKS6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:50-1:100

Target information (Q9UKS6)


Gene Name	PACSIN3
About Gene Symbol	PACSIN3 (Protein Kinase C And Casein Kinase Substrate In Neurons 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PACSIN3 in HeLa, A549 lysates using PACSIN3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using PACSIN3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PACSIN3
(green) in Hela using PACSIN3 antibody, and
DAPI(blue)