


ING1 Rabbit mAb

ABZ1764


Product Overview

Name	ING1 Rabbit mAb
Accession(Primary)	Q9UK53
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (Q9UK53)

Gene Name ING1**About Gene Symbol**

Western blot analysis of ING1 in C6, Raw264.7 lysates using ING1 antibody.



Western blot analysis of ING1 in HeLa, A549, HL-60, U251, U87-MG lysates using ING1 antibody.