

LEF1 Rabbit mAb


ABZ1766

Product Overview



Name	LEF1 Rabbit mAb
Accession(Primary)	Q9UJU2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 25-58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:50-1:200; FC-1:20-1:50

Target information (Q9UJU2)

Gene Name	LEF1
About Gene Symbol	



Immunocytochemistry analysis of LEF1 (green) in Jurkat using LEF1 antibody, and DAPI (blue).



Western blot analysis of LEF1 in Jurkat lysates using LEF1 antibody.