

Sall4 Rabbit mAb

ABZ1767

Product Overview

Name	Sall4 Rabbit mAb
Accession(Primary)	Q9UJQ4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 112 kDa; Observed MW: 112 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (Q9UJQ4)

Gene Name SALL4**About Gene Symbol**

