

DAPK2 Rabbit mAb


ABZ1768

Product Overview

Name	DAPK2 Rabbit mAb
Accession(Primary)	Q9UIK4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:200-1:1000

Target information (Q9UIK4)

Gene Name	DAPK2
About Gene Symbol	



Western blot analysis of DAPK2 in rat spleen
lysates using DAPK2 antibody.