

## SUN2 Rabbit mAb


ABZ1769

### Product Overview

|                             |                                                                            |
|-----------------------------|----------------------------------------------------------------------------|
| <b>Name</b>                 | SUN2 Rabbit mAb                                                            |
| <b>Accession(Primary)</b>   | Q9UH99                                                                     |
| <b>Host Species</b>         | Rabbit                                                                     |
| <b>Clonality</b>            | Monoclonal                                                                 |
| <b>Isotype</b>              | IgG                                                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                          |
| <b>Calculated MW</b>        | Calculated MW: 80 kDa; Observed MW: 80 kDa                                 |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:1000;<br>ICC/IF-1:200-1:2000; FC-1:20-1:50 |


### Target information (Q9UH99)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | SUN2 |
| <b>About Gene Symbol</b> |      |



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of SUN2 in HepG2 and Jurkat lysates using SUN2 antibody.



Western blot analysis of SUN2 in Jurkat lysates using SUN2 antibody.