

SPAK Rabbit mAb


ABZ1770

Product Overview


Name	SPAK Rabbit mAb
Accession(Primary)	Q9UEW8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 59 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q9UEW8)

Gene Name	STK39
About Gene Symbol	STK39 (Serine/Threonine Kinase 39) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SPAK in HeLa, A549, HL-60, U2OS lysates using SPAK antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using SPAK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.