


Wnt16 Rabbit mAb

ABZ1772


Product Overview

Name	Wnt16 Rabbit mAb
Accession(Primary)	Q9UBV4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9UBV4)

Gene Name WNT16**About Gene Symbol**

Western blot analysis of Wnt16 in mouse brain
lysates using Wnt16 antibody.



Western blot analysis of Wnt16 in K562, rat Brain lysates using Wnt16 antibody.