

NXF1 Rabbit mAb


ABZ1773

Product Overview



Name	NXF1 Rabbit mAb
Accession(Primary)	Q9UBU9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:10-1:100; IHC-1:50-1:100

Target information (Q9UBU9)


Gene Name	NXF1
About Gene Symbol	



Immunocytochemistry analysis of NXF1 (green) in Jurkat using NXF1 antibody, and DAPI (blue).



Western blot analysis of NXF1 in K562 lysates using NXF1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using NXF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.