

UBA2 Rabbit mAb


ABZ1776

Product Overview


Name	UBA2 Rabbit mAb
Accession(Primary)	Q9UBT2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (Q9UBT2)

Gene Name	UBA2
About Gene Symbol	



Immunocytochemistry analysis of UBA2 (green) in 293 using UBA2 antibody, and DAPI (blue).



Western blot analysis of SAE2 / UBA2 in rat Brain, C6, HeLa lysates using SAE2 / UBA2 antibody.