


DcR2 Rabbit mAb

ABZ1779


Product Overview

Name	DcR2 Rabbit mAb
Accession(Primary)	Q9UBN6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000

Target information (Q9UBN6)

Gene Name TNFRSF10D**About Gene Symbol**

Western blot analysis of DcR2 in Jurkat lysates using DcR2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DcR2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.