

PGC1 alpha Rabbit mAb

ABZ1780


Product Overview

Name	PGC1 alpha Rabbit mAb
Accession(Primary)	Q9UBK2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 113 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9UBK2)

Gene Name PPARGC1A

About Gene Symbol



Western blot analysis of PGC1 alpha in HeLa, NIH 3T3, and C6 lysates using PGC1 alpha antibody.

