

MAGEC2 Rabbit mAb


ABZ1781

Product Overview


Name	MAGEC2 Rabbit mAb
Accession(Primary)	Q9UBF1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q9UBF1)

Gene Name	MAGEC2
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse testis using MAGEC2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of MAGEE1 in K562 lysates using MAGEE1 antibody.