


SAE1 Rabbit mAb

ABZ1782


Product Overview

Name	SAE1 Rabbit mAb
Accession(Primary)	Q9UBE0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100

Target information (Q9UBE0)

Gene Name SAE1**About Gene Symbol**

Immunocytochemistry analysis of SAE1 (green) in Hela using SAE1 antibody, and DAPI(blue).



Western blot analysis of SAE1 in HeLa, Jurkat
lysates using SAE1 antibody.