



EPS15R Rabbit mAb

ABZ1783


Product Overview

Name	EPS15R Rabbit mAb
Accession(Primary)	Q9UBC2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 94 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:2000-1:5000; ICC/IF-1:50-1:100

Target information (Q9UBC2)

Gene Name EPS15L1**About Gene Symbol**

Immunocytochemistry analysis of EPS15R (green)
in HeLa using EPS15R antibody, and DAPI(blue).



Western blot analysis of EPS15R in A549, HL-60
lysates using EPS15R antibody.