

TLR2 Rabbit mAb


ABZ1785

Product Overview


Name	TLR2 Rabbit mAb
Accession(Primary)	Q9QUN7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 89 kDa; Observed MW: 89 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (Q9QUN7)


Gene Name	Tlr2
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse lung using Toll Like Receptor 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TLR2 in Raw264.7 lysates using TLR2 antibody.



Immunocytochemistry analysis of TLR2 (green) in Raw264.7 using TLR2 antibody, and DAPI (blue)