

MST4 Rabbit mAb


ABZ1786

Product Overview


Name	MST4 Rabbit mAb
Accession(Primary)	Q9P289
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500; FC-1:100-1:1000; IP-1:20-1:50

Target information (Q9P289)

Gene Name	STK26
About Gene Symbol	



Western blot analysis of MST4 in K562, C6 lysates using MST4 antibody.



Immunocytochemistry analysis of MST4(green) in Hela using MST4 antibody, and DAPI(blue)