

PIM2 Rabbit mAb


ABZ1787

Product Overview

Name	PIM2 Rabbit mAb
Accession(Primary)	Q9P1W9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 40,38,34 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9P1W9)

Gene Name	PIM2
About Gene Symbol	



Western blot analysis of PIM2 in A549, HL-60, U251 lysates using PIM2 antibody.