

OR1A1 Rabbit mAb

ABZ1788


Product Overview


Name	OR1A1 Rabbit mAb
Accession(Primary)	Q9P1Q5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20

Target information (Q9P1Q5)

Gene Name OR1A1

About Gene Symbol





Western blot analysis of OR1A1 in CHO-K1, C6, Jurkat lysates using OR1A1 antibody.