

NDUFAF4 Rabbit mAb


ABZ1789

Product Overview


Name	NDUFAF4 Rabbit mAb
Accession(Primary)	Q9P032
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; IP-1:50-1:100

Target information (Q9P032)

Gene Name	NDUFAF4
About Gene Symbol	



Immunocytochemistry analysis of NDUFAF4 (green) in HCT116 using NDUFAF4 antibody, and DAPI (blue).



Western blot analysis of NDUFAF4 in K562, HeLa lysates using NDUFAF4 antibody.