

CDT2 Rabbit mAb


ABZ1791

Product Overview


Name	CDT2 Rabbit mAb
Accession(Primary)	Q9NZJ0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:100-1:200

Target information (Q9NZJ0)

Gene Name	DTL
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using CDT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CDT2 in HeLa lysates using CDT2 antibody.