

PEN2 Rabbit mAb


ABZ1794

Product Overview


Name	PEN2 Rabbit mAb
Accession(Primary)	Q9NZ42
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:100-1:500; ICC/IF-1:50-1:100

Target information (Q9NZ42)


Gene Name	PSENN
About Gene Symbol	




Immunocytochemistry analysis of PEN2 (green) in Hela using PEN2 antibody, and DAPI (blue).



Western blot analysis of PEN2 in Raw264.7, mouse brain lysates using PEN2 antibody.



Western blot analysis of PEN2 in HeLa, A549, HL-60, C6 lysates using PEN2 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PEN2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.