

ETAA1 Rabbit mAb

ABZ1795


Product Overview

Name	ETAA1 Rabbit mAb
Accession(Primary)	Q9NY74
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 103 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9NY74)

Gene Name ETAA1

About Gene Symbol



Western blot analysis of ETAA1 in K562 lysates using ETAA1 antibody.

