

## ALKBH4 Rabbit mAb

ABZ1796


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | ALKBH4 Rabbit mAb                          |
| <b>Accession(Primary)</b>   | Q9NXW9                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 34 kDa; Observed MW: 34 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20-1:50              |

### Target information (Q9NXW9)

**Gene Name** ALKBH4

**About Gene Symbol**



Western blot analysis of ALKBH4 in K562, HeLa lysates using ALKBH4 antibody.

