

HIF Prolyl Hydroxylases Rabbit mAb


ABZ1797

Product Overview


Name	HIF Prolyl Hydroxylases Rabbit mAb
Accession(Primary)	Q9NXG6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (Q9NXG6)

Gene Name	P4HTM
About Gene Symbol	



Western blot analysis of HIF Prolyl Hydroxylases in HeLa, A549, HL-60, U251, U87-MG lysates using HIF Prolyl Hydroxylases antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HIF Prolyl Hydroxylases antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.