

IRAK4 Rabbit mAb


ABZ1799

Product Overview


Name	IRAK4 Rabbit mAb
Accession(Primary)	Q9NWZ3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:50-1:100; IHC-1:50-1:100

Target information (Q9NWZ3)


Gene Name	IRAK4
About Gene Symbol	



Western blot analysis of IRAK4 in Jurkat, HeLa lysates using IRAK4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IRAK4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of IRAK4 (green)
in hela using IRAK4 antibody, and DAPI(blue)