

PRMT7 Rabbit mAb


ABZ1801

Product Overview


Name	PRMT7 Rabbit mAb
Accession(Primary)	Q9NVM4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9NVM4)

Gene Name	PRMT7
About Gene Symbol	PRMT7 (Protein Arginine Methyltransferase 7) is a methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRMT7 in rat muscle, rat pancreas lysates using PRMT7 antibody.



Western blot analysis of PRMT7 in K562, 3T3, HeLa lysates using PRMT7 antibody.