


POT1 Rabbit mAb

ABZ1802

Product Overview

Name	POT1 Rabbit mAb
Accession(Primary)	Q9NUX5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9NUX5)

Gene Name POT1**About Gene Symbol**

Western blot analysis of POT1 in K562, C6, 3T3 lysates using POT1 antibody.

