

SIRT3 Rabbit mAb


ABZ1803

Product Overview


Name	SIRT3 Rabbit mAb
Accession(Primary)	Q9NTG7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:2000

Target information (Q9NTG7)

Gene Name	SIRT3
About Gene Symbol	SIRT3 (Sirtuin 3) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets



Western blot analysis of SIRT3 in Jurkat, HeLa lysates using SIRT3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using SIRT3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.