

TREX1 Rabbit mAb


ABZ1805

Product Overview


Name	TREX1 Rabbit mAb
Accession(Primary)	Q9NSU2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (Q9NSU2)


Gene Name	TREX1
About Gene Symbol	TREX1 (Three Prime Repair Exonuclease 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of TREX1 (green) in HeLa using TREX1 antibody, and DAPI (blue).



Western blot analysis of TREX1 in HeLa lysates using TREX1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using TREX1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.