

PLSCR3 Rabbit mAb


ABZ1806

Product Overview


Name	PLSCR3 Rabbit mAb
Accession(Primary)	Q9NRY6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9NRY6)

Gene Name	PLSCR3
About Gene Symbol	PLSCR3 (Phospholipid Scramblase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PLSCR3 in rat lung, mouse lung lysates using PLSCR3 antibody.



Western blot analysis of PLSCR3 in K562, HeLa lysates using PLSCR3 antibody.