

Striatin 4 Rabbit mAb


ABZ1807

Product Overview

Name	Striatin 4 Rabbit mAb
Accession(Primary)	Q9NRL3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:10-1:100; IP-1:10-1:100

Target information (Q9NRL3)

Gene Name	STRN4
About Gene Symbol	STRN4 (Striatin 4) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Striatin 4 in HeLa lysates using Striatin 4 antibody.