


MMP26 Rabbit mAb

ABZ1809

Product Overview

Name	MMP26 Rabbit mAb
Accession(Primary)	Q9NRE1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9NRE1)

Gene Name MMP26**About Gene Symbol**

Western blot analysis of MMP26 in Jurkat, C6, CHO-K1 lysates using MMP26 antibody.

